



Chemical Compatibility of Eppendorf Consumables

Clemson University Safety Retreat
August 11th, 2022

Eddie Bondo
Senior Territory Manager
Western NC, all of SC

1



Well-Known Eppendorf Products



2

eppendorf

Chemical Stability Guide - Consumables

eppendorf

USERGUIDE No. 23


Chemical Stability of Consumables

Kerstin Bock, Eppendorf AG, Hamburg, Germany

Abstract
Consumables such as Eppendorf Tubes® and pipette tips have to withstand the different kinds of stress and strain encountered in every day laboratory work, e.g. when centrifuging and autoclaving, and at the same time are exposed to chemical reactions. The level of stability to these effects is, to a large degree, determined by their material and its processing during manufacture. As a leading manufacturer of labware products, Eppendorf uses PP and PE granulates specifically suited for the laboratory applications and manufacturing process of its micro test tubes, Concomet® and tips.

Introduction
Chemical interactions with consumables
Chemicals interact with consumables in a number of different ways. These include:
1. The chemical reaction with the polymer, e.g. by oxidation, mediated with the functional groups, and by depolymerisation.
2. Changing the specific chemical surface to cause leaching the surface into the solution, material to induce swelling and a change in the material's mechanical properties.
3. The formation of inclusion cracks resulting from mechanical stress and surface active substances (detergents, salts). The stress can be caused by repeated freeze-thaw cycles or as a consequence of an external process (autoclaving).
Potential applications
In many applications it is necessary to dispense acids, bases, solvents and salt solutions with varying densities, vapor pressures and viscosities, whereby their properties differ considerably from that of water. Such substances influence the chemical resistance of the materials with which contact is made. After using a dispensing system, two questions must thus be answered:
1. How great are the dispensing errors resulting from the physical properties of the liquid?
2. Is the pipette tip or the Concomet® during dispensing despite the chemical resistance and the pipette or dispenser resistant for a specified period of time?

Chemical stability
Consumables are generally made of polypropylene (PP) or polyethylene (PE). For both thermoplastics the relevant literature contains extensive information about chemical resistance, in terms of a high level of chemical resistance for both PP and PE (1). According to this reference work, chemical stability is not attributed to plastic, for example, even though it may be designed to withstand contact with a chemical for a period of weeks. An example is this case for the storage of acetone in plastic tubes; it does not affect conditions present for the most three contact times combined with dispensing. Liquids to which PP or PE are only conditionally resistant according to the literature such as inorganic acids and various organic solvents such as acetone, diethyl ether, chloroform or toluene can be dispensed without a problem. The present research work is that the contact time is increased, i.e. the plastic is to be used for unlimited contact periods on the basis of an additional test.
As regards the resistance of pipette and dispenser tips, it is a basic rule that the expected dispensing error will increase the longer a consumable is used for the above



3

eppendorf

Using this Guide...

Materials in selected Eppendorf products

	PE Polyethylene	PP Polypropylene	PC Polycarbonat	autoclavable
epT.I.P.S.® All tips	■	■	■	■
epT.I.P.S.® Box, epT.I.P.S.® Reloads		■		■
Box: min. autoclavable at least 100 times		■		■
Reloads: autoclavable		■		■
ep Dualfilter T.I.P.S.® All tips	■	■		■
Filter: hydrophobic	■			■
Mastertip®		■		■

Eppendorf Tubes® (3810X, Safe Lock, Eppendorf LoBind®)
All tubes
When lid open

Eppendorf Plates® (DWP, MTP, LoBind)
If unsealed, but the dimensional stability can be affected

Combitips advanced®
Cylinder
Piston
Adapter for 25 mL and 50 mL: Polybutylene terephthalate (PBT) min. autoclavable at least 100 times

4

Using this Guide...

Chemical resistance of PP and PE

Chemicals	Concentration in %	Other names	Temperature °C	PE ¹	PP ²	PC ³	Steam pressure at 20 °C (hPa)	Density at 20 °C (mg/μL)	Viscosity at 20 °C (mPas)	Remarks
Acetaldehyde	40	Ethanal	20 60	0 3	1 1	3	1006	0.78		
Acetic acid	25-60	Ethane acid	40 60	1 1	1 1	0		1.06	1.22	
Acetone	100	Dimethylketone	20 60	2 3	1 2	3	246	0.79		
Ammonium hydroxide	30		60	1	2	0	483	0.89		

1 = resistant; the material does not change even after longer contact with the substance.
 2 = conditionally resistant; if the material gets in contact with the substance only for a short time it does not change.
 3 = non-resistant; the material already changes after a short contact with the substance.
 0 = no existing value.

5

Thanks!

For questions, please come visit me at my table...

6